

Station A: MagMAX Viral/Pathogen Nucleic Acid Isolation


Code parameters:

- Change the sample number on line 14 (default is 8, max is 96)
- Change the sample volume (µl) on line 15 (default is 200)
- Tip rack tracking can be changed from False to True on line 16 (default is False)

Pipettes:

- P1000 single channel on the right mount
- P20 multi channel on the left mount

Deck Layout:



Labware and module requirements:

- ≤ 96 x 1.5 - 2mL tubes [input - samples]
- 4 x 24 Tube rack [holds 1.5 - 2mL tubes with samples]

- 3 x 1000µl filter tip racks
- 1 x 6 Tube rack [**holds 1 50mL tube with Binding Bead Mix in A1**]
- 1 x 50mL tube [**Binding Bead Mix**]
- 1 x 2mL deep well plate [**output - destination plate**]
- 1 x 20µl filter tip racks
- 1 x 1.5-2mL tube [**for creating the internal control & proteinase k mix**]
- 1 x 96 well aluminum block [**holds 1 PCR strip tube in column 1**]
- 1 x 200µl PCR strip tubes [**holds internal control & Proteinase K mix**]

Volume requirements:

Note: the below volumes account for a dead volume - the dead volume can be adjusted depending on the calibration of the pipette to the labware, but it's best to have a dead volume of at least 10%


Reagent	Volume per sample	Volume for 8 samples	Volume for 48 samples	Volume for 96 samples
Binding Bead Mix	275µl	3,300µl	14.3mL	27.5mL
Proteinase K & Internal Control mix	10µl	120µl	520µl	1000µl

Before you begin:

1. Pre-cool the Temperature Module in the Opentrons App to 4°C
2. Load the samples onto the tube racks starting with filling up slot 2, slot 3, slot 5, slot 6.
3. Create the **Binding Bead Mix for slot 7**

Reagent	Volume per sample	Volume for 8 samples	Volume for 48 samples	Volume for 96 samples
Binding Solution	265µl	3,180µl	13.7mL	26.5mL
Total Nucleic Acid Magnetic Beads	10µl	120µl	520µl	1mL


4. Add the binding bead mix in the 50mL tube to the rack in **Slot 7**



5. Create the **Proteinase K and Internal Control mix** in a 1.5-2ml tube:

Reagent	Volume per sample	Volume for 8 samples	Volume for 48 samples	Volume for 96 samples
Proteinase K	5µl	60µl	260µl	500µl
Internal Control	5µl	60µl	260µl	500µl

6. Divide the total volume of the Proteinase K and Internal Control mix by 8 (example: 96 samples would have 125µl in each PCR tube) and add it to a 200µl PCR strip tube to be placed in **Slot 1** on top of a 96 well aluminum block in column 1



The final destination plate of lysates will be in slot 4 in the deep well plate. Once the run is finished, move the deep well plate to Station B to complete the remainder of the extraction protocol.